

réjuvapen NXT[®]

**TRANSFORM YOUR
MICRONEEDLING**

rejuvapennxt.com

WHAT IS MICRONEEDLING?

Microneedling is a widely-practiced, minimally invasive way to rejuvenate skin using the body's natural ability to heal itself. Fine sterile needles create small punctures in the epidermal and dermal layers of the skin.

This leads to the release of growth factors and stimulates the formation of new collagen, elastin, and neovascularization in the dermis.

Microneedling results in healthier, younger-looking skin.



WHY REJUVAPEN NXT?



Rejuvapen NXT® is an FDA cleared Class II device to treat facial wrinkles. With its advanced design, the Rejuvapen NXT® delivers both improved profitability and exceptional results.

The system advances microneedling with its **nine speed settings, precise Infinite Dial, and brushless digital motor.**

With Rejuvapen NXT®, providers now have more flexibility to customize each patient's treatment to produce **consistent clinical results.**




Tru Vertical® Microneedling

Tru Vertical® microneedling results in more effective healing while **keeping the overall structure of the skin intact.**

The Tru Vertical® technique uses **a gentle stamping-like motion instead of side to side or circular movement.** This simple light touch against the skin is repeated until the area of the skin has been fully covered.

The Tru Vertical® method is taught in the **live online microneedling instruction included with the device.**




PRODUCT SPECS



PATENTED ISO SEAL

PRECISE INFINITE DIAL

**LIGHTWEIGHT
MEDICAL GRADE
ALUMINUM BODY**





Infinite Dial

Precision is everything: Depth settings of 0.1 to 2.5mm for more precise treatments.

Smart Touch

Nine speed settings to better
customize treatments



Fluid barrier between tip
and device

12 stainless steel tips (33-
gauge) for optimal results

More secure threaded attachment

ISO TipTM

U.S. Patent . No. 10,945,759

Scalloped edge for improved safety
with no suction or clogging



Fluid-resistant barrier
between tip and device

12 stainless steel tips (33-
gauge) for optimal results

4x The treatment area!



EXCEL

PATIENT BENEFITS

FASTER RECOVERY

True vertical punctures result in less trauma to surrounding skin

AFFORDABLE

With the lowest per treatment cost, you have more room to offer specials or pass savings to your patients

LESS DISCOMFORT

Personalized speed and depth settings ensure comfort





RESULTS MATTER

Clinically proven to be safe and effective for all skin types

Optimized for stimulating collagen and elastin production

Improves the appearance of fine lines and wrinkles

IMAGES FROM OUR PROVIDERS



Photos courtesy of Cassie Bosco | Amazing Skin, LLC | Allentown, PA



Photos courtesy of Christl Colven | Siempre SkinCare | Boise, ID



Photos courtesy of Mary Richardson | Permanent Makeup by Mary | Tampa, FL



PROVIDER BENEFITS


*OF FDA CEARED MICRONEEDLING DEVICES

Lowest per treatment cost* and highest ROI.

Proven safe and effective in more than 1,000,000 treatments.

Live online training included at no charge Better results and higher patient satisfaction.

Plug-in system for full power and consistent output.



”

LAURA WOODS FNP-C*Lafayette, Louisiana*

Rejuvapen NXT has been a vital tool in my practice to accelerate patients results when paired with medical grade skincare. The Tru-Vertical technique has provided me more control of the treatment, and less downtime for my patients with amazing results!”



”

ROBERT J. ALLEN. M. D.*Arizona Wellness Center*

"We love using the Rejuvapen NXT for microneedling! Its ability to be dialed to precisely the depth we need, as well as the option to adjust speed at which the needles penetrate the skin, make it easy to use. We've seen great results, especially in the reduction of fine lines and wrinkles. Our patients are thrilled to see their skin improve with each treatment."



SOCIAL MEDIA & MARKETING SUPPORT

30 DAY ADVERTISING AND
MARKETING PLAN

PRE-MADE CANVA TEMPLATES FOR
SOCIAL MEDIA AND WEB

PATIENT BROCHURES FOR IN
PRACTICE MARKETING





THANK YOU!

CONTACT US



Phone

866-590-5533



Mail

sales@refineusa.com



Website

www.rejuvapennxt.com



Address

340 3rd Ave S C, Jacksonville Beach, FL 32250